

Stryker

800 & 9000 Series Autopsy Saws

810



Armature - CAST 044 (See Armature/Field Assembly)
Armature Bearing Front - CAST026
Armature Bearing Rear - CAST027
Armature/Field Assembly - CAST133
Armature Housing Spring - CAST064
Blade Adapter Autopsy - CAST116
Blade Sectioned Autopsy - CAST115
Blade Sectioned Autopsy with Shaft - CAST121
Brush - CAST023
Brush Cover (800 Series) CAST022
CAM - CAST003
CAM Bearing - CAST025
Field 115 Volt - CAST041 (See Armature/Field Assembly)

9004-210



Internal Retaining Ring - CAST074
Link Compatible with pin Drive Shaft - CAST008
Oil Seal - CAST029
On/Off Switch - CAST032
Power Cord - CAST048
Set Screw - CAST134
Strain Relief Fitting - CAST007
Strain Relief Spring - CAST070

8381-210

